

TECHNICAL DATA SHEET

T HEXX DRY-E TEAT SEALANT

Ready to use drying off flexible barrier for use as part of the drying off routine either in addition to or instead of traditional dry cow therapy or internal teat sealants.

PRODUCT BENEFITS:

- Ideal complement to traditional dry cow therapy.
- T-Hexx® Dry-E assists the natural keratin plug which is the teat's defence mechanism.
- Forms a physical flexible barrier on and in the teat for up to 7 days.
- Quick drying within 5-10 minutes and is resistant to water, manure & urine.

USAGE GUIDE:

- T Hexx Dry-E should be applied at drying off and again 7 days prior to calving if required.
- This product is also ideal for use on in-calf heifers, especially in areas with high fly populations.

INSTRUCTIONS FOR USE:

Postmilking use:

1. Fully milk out quarter.
2. Disinfect and dry teats and infuse drying off tubes/internal teat sealants.
3. Dip each teat using T-HEXX DRY-E sealant.
4. Hold cup under teat to catch excess drips for a few seconds.
5. Wait for 5 minutes before releasing animal until material is mostly dry.
6. Do not pour unused product back into container.

Shelf life:

1 year from date of manufacture.

Storage:

Store between 0°C to 25°C.

Keep out of direct sunlight.

For health and safety information, please refer to the Safety Data Sheet.



SIZES AVAILABLE:

- 950ml (95 cows) - Code P415